TGN-RSG70

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CH-REP. CH-Importeur: INNOVAhelp GmbH, Thunstr. 5a, CH-3110 Münsingen

Dear Customer,

We thank you for your confidence in our company by purchasing one of our products. We are convinced that our trapeze handle entirely meets the requirements of the daily nursing routines. Our trapeze handle was developed according to the highest medical standards, using the most advanced technologies and with the thoroughness and the team spirit of our employees. We are confident that our trapeze handle makes everyday life of nursing staff and people in need of care easier and more pleasant. Please take note of instructions listed below and consider specially the safety instructions to guarantee the proper use of our product.

01. Purpose and Indication

These operating instructions are intended to qualified nursing staffs but also to any person who is in charge of the initial and further operation of this product. All drawings and graphics are not drawn to scale and may deviate significantly. The product is intended to be used only in the field of nursing, hospital and rehabilitation. It serves as an accessory for nursing/hospital beds compensating limitations imposed by disabilities and making nursing tasks. Patients can stand up with the help of the handle independently, provided that they have the strength to do so. The height of the handle can be quickly and easily adapted individually to the patient.

Other applications must be approved by the manufacturer in advance. The product is to be used as a working tool for nursing staff and is subject to the regulations of the employer's Liability Insurance Association (in German: Berufsgenossenschaft BG).

02. Contraindication


The handle is not suitable for all applications that differ from the intended use and or may lead to danger. These include, for example, exceeding the permitted max. load (see: *Safety Instructions and Remaining Risks* and *Technical Data*) or the operation of the handle by untrained staff. Please note that the handle is not suitable as a sport equipment for rehabilitation exercises.

03. Safety Instructions and Remaining Risks

- These operating instructions are to be read prior to use of the handle.
- The handling of this product requires exact knowledge and observation of the instructions listed below to avoid operating mistakes and grant a failure-free operation. Keep this instruction manual for the operating personnel always accessible until the disposal of the product.
- Do not use the handle against the intended purpose and avoid improper use.
- Before using the handle the user needs to familiarize himself with the functional reliability of the handle and verify the proper condition of the trapeze bar.
- Check the product for damages regularly (see: *Maintenance and Inspection*).
- Replace the handle immediately if any damage to the handle and or its components are visible. Observe the cleaning and disinfection specifications.
- Pay attention to the correct adjusted height of the handle in relation to the patient.
- For safety reasons, avoid the permeation of small foreign bodies and/or liquids in the handle (see: *Care, Cleaning and Disinfecting*).
- Keep the handle away from strong heat sources such as cigarettes or open flames and observe the specified ambient conditions below.
- The handle must not be often rotated around the vertical axis (over 180°) in suspended state.
- Do not exceed the safe working load of 75 kg (exception: TGN-RSGP70S 90 kg).
- If serious incidents occur in connection with this product, the incidents must be reported to Rolko Kohlgrüber GmbH and to the competent authority of the EU Member State in which the user or patient is based.

04. Operation

Read this instructions manual carefully and check if the product is complete or has any damages. The length of the belt may need to be preset before first use. An initial cleaning of the handle is necessary, in compliance with the permitted cleaning substances and disinfectants (see: *Care, Cleaning and Disinfecting*), this also applies for the further uses of the product. Perform the mounting on the trapeze bar as described in the following instructions.

• Attachment to the trapeze bar


Hang the handle with the loop on the trapeze bar. The slip guard has to properly between belt strap and trapeze bar.

The attachment is only allowed on a trapeze bar, which is intended for this handle. Please also consider the instructions notes of the trapeze bar.

• Height adjustment

The height adjustment has to be carried out by adjusting the belt strap with the plastic buckle. The handle must always be relieved.

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Do not operate the handle under load! The handle must always be relieved.

• Application and use

Set the ideal position of the handle for the patient so that he or she can use it as an aid for sitting up or moving in bed. The ergonomically shaped grip surface has an anti-slip structure. The handle may be used with one or both hands.

Do not overload the handle max. 75 kg (exception: TGN-RSGP70S 90 kg)!

• Park position when not in use

While the handle is not being used, you may hang it on the trapeze bar. It must be ensured that the handle can not accidentally slip off the trapeze bar.

• Functional check

Check the function of the product according to section 6: *Function description - Height adjustment.*

05. Care, Cleaning and Disinfecting

Before cleaning and disinfecting the handle, remove it from the bed gallows / trapeze bar. Please consider the following guidelines when cleaning.

The cleaning and disinfection of the product must not be carried out in pressure chambers / autoclaves, sterilisation, cleaning or disinfection machines!

Clean and disinfect the handle before the first use and also before each reuse. Disinfect the handle on a regular basis and at least prior to each change of patient. The shape of the handle is designed to facilitate the cleaning. In addition, the handle is made of polypropylene and thus it has a high chemical resistance, which makes the hygiene process easier. There were no additional biocidal agents used during the manufacturing of the product so that allergy and resistance risks were not increased.

• Care and Cleaning

For an easy cleaning the trapeze handle may be wiped clean with a damp cloth and if necessary with the addition of a mild detergent. Use lukewarm water and if necessary, non-abrasive, pH-neutralized and commercial household cleaners. Please also note the following section *Disinfection / Sterilisation*. After each cleaning and disinfection please let it dry completely so that no residual moisture remains in the webbing or other parts.

• Disinfection / Sterilisation

We recommend the detergents and disinfectants listed below to clean the handle, after consultation with the manufacturers. The usage must be carried out only according to manufacturer specifications! Observe the corresponding product information, in particular regarding the dosage and a possible hazard in combination with other cleanings / disinfection agents. Before using other detergents or disinfectants please consult the agent manufacturers in any case. Due to the variety of detergents and disinfectants available on the market, as well as constant innovation, we can not guarantee any kind of compatibility. In case of doubt please contact the manufacturer, since a misapplication may lead to a premature wear and failure of the material presenting a danger to the user.

- From the field of formaldehydes and/or other aldehydes respective derivatives: Aldasan 2000, B5, Budenat Acute D 441, Desifor Protect, Desomed Perfekt, hygienicVIRUZID, Kohrsolin extra, Lysoform, Lysoformin, Melsitt, Minutil, Nüscosept
- From the field of percompounds: APESIN AP100, APESIN AP 100 PLUS, DISMOZON PLUS, DISMOZON PUR, INCIDIN ACTIVE, PERFORM, TERRALIN PAA, WOFASTERIL
- From the field of carboic acid (derivate): AMOCID

The following substances may **not** be used, since they alter the material properties, leading to premature wear, which may cause a potential material failure and thus hazards to the user:


- All phenol or phenol derivatives not mentioned in these instructions for use and other aromatic hydrocarbons.
- Chlorine organ. or inorganic substances with active chlorine.
- Alkalis.
- All formaldehyde and / or other aldehydes or derivatives not mentioned in these instructions for use.

For safety reasons, the permeation of small foreign bodies or liquids (not allowed, e.g. solvents, oils, alkalines, acidic media, cleaners etc.) in any part of the handle must be avoided.

06. Maintenance and Servicing

The lifetime of the handle including the belt strap is influenced by many factors, such as temperature, humidity, UV radiation, exposure to chemicals, mechanical stress, etc. Also, a careless handling may cause wearout and premature aging of the plastic parts, a manufacturer has no influence in this process. Therefore, we recommend to exchange the product after every 5 years as a precaution, even if no signs of external damage exist. The recommended service life, as well as the maintenance intervals stated below are valid from the date of commissioning.

Manufacturer

Follow the instructions

Production date

Item no.

Batch

Medical device

Unique device identification

- **Maintenance intervals**

| | |
|-----------------------------------|-----------|
| Functional test: | Quarterly |
| Control of plastic parts: | Quarterly |
| Control of webbing and stitching: | Quarterly |
| Control of anti-slip protection: | Annually |

If any kind of defects or problems appears during functional tests or maintenance, the product must be replaced immediately.

- **General test**

Inspect the trapeze handle regularly regarding safe functionality and visual impairments. All visible parts must be inspected regularly for wearout or damage marks, for example, by scratches, discolorations, cracks etc. The webbing and stitching must not be cut, nor look twisted or frayed.

If the trapeze handle and/or its components are damaged or show any signs of malfunction, the product must be replaced immediately!

07. Storage and Reutilization

Before storing the handle the following steps need to be performed:

1. Conduct a function test before storing the product.
2. Cover or wrap the handle so that it does not get dirt or damaged. Protect from prolonged/permanent UV light.
3. Due to the maintenance periods attach a date of storage to the product.
4. Ensure that this operating instructions are present.
5. Store the handle only in closed rooms with a temperature between +10 and +40 °C and with about max. humidity of 80%. Before reusing it please conduct the cleaning and disinfection steps as described under section 7: *Care, Cleaning and Disinfecting*.

08. Disposal

All parts of the handle must be disposed of safely with household waste. To do this, disassemble the grab handles into its individual parts to enable the materials to be sorted.

09. Technical Specifications

Dimensions of the handle (height x width): ca. 155 x 300 mm

Height adjustment range (from handle grip area to the upper edge of the trapeze bar): ca. 320 bis 500 mm (may vary depending on the model)

Max. load capacity: 75 kg (exception: TGN-RSGP70S 90 kg)

Own weight: ca. 286 g (may vary depending on the model)

Materials used: PP, POM, PES and PVC

Permitted ambient conditions:

Temperature between +10 °C bis +40 °C

Relative humidity up to max. 80%

As an accessory for medical products of class 1, this product has a CE mark in accordance with regulation (EU) 2017/745 (MDR) and thus complies with all German and European safety requirements, that apply to the product. The product fulfils the corresponding sections of standard specification DIN EN 60601-2-52. The EU declaration of conformity is available for download at instructions.rolko.com.

Operating Instructions: Medical Trapeze Handle with Belt Strap: TGN-RSG70, TGN-RSG70-V, TGN-RSG100-SO, TGN-RSG250, TGN-RSG90, TGN-RSGP70S / 2024-06-27 / Rev. 12 / Basis-UDI-DI: 426067605TGN0000000000UN

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Reprint, also in extracts, only with prior permission by the editor. All rights reserved.

The delivered product may differ from the illustrations in this manual.

All technical data are characteristic and subject to construction and manufacturing conditioned tolerances.

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© Rolko Kohlgrüber GmbH, Industriestr. 14 - 16, 33829 Borgholzhausen